



Model	Pharm260	Pharm360	Pharm610	Pharm890	Pharm1378
Capacity(L)	260	360	610	890	1378
Voltage/Frequency	220/50	220/50	220/50	220/50	220/50
Power(W)	340	430	490	900	800
Temp. range(°C)	2~8	2~8	2~8	2~8	2~8
Door	Single-layer glass	Single-layer glass	Single-layer glass	Double-layer glass	Double-layer glass
Shelves/Inner door	4/-	5/-	6/-	12/-	12/-
Internal DI(mm)(W*D*H)	550*460*1065	550*460*1340	680*620*1400	1030*590*1425	1320*700*1500
Exterior DI(mm)(W*D*H)	620*655*1720	620*655*1955	780*840*1960	1125*740*1980	1440*875*1980
Packing DI (mm)(W*D*H)	700*715*1810	700*715*2080	875*937*2097	1260*850*2150	1600*1010*2100
Net/Gross Weight(kg)	88/100	99/114	204/227	227/279	310/360

Model	Sun60
Storage Capacity(Litres)	21
Gross Capacity(Litres)	60
Ambient Temperature	10C~40C
Temp. Range (C)	2~8
Autonomy Time at 30C	130 h
Hold Over Time During Power Cut	110-120 h
Internal DI(mm)(W*D*H)	240*445*610
Exterior DI(mm)(W*D*H)	788*745*875
Net/Gross Weight(kg)	78/98



Model	BB158	BB258	BB358	BB608	BB890	BB1308	BB1308K
Capacity(L)	158	258	358	608	890	1308	1308K
VAC /HZ	220/50	220/50	220/50	220/50	220/50	220/50	220/50
Power(W)	350	430	460	490	900	850	850
(°C)	4±1	4±1	4±1	4±1	4±1.5	4±1.5	4±1.5
Temp.control	Microprocessor control				Microprocessor control		
Ambient Temp.	10°C~38°C	10°C~38°C	10°C~38°C	10°C~38°C	10°C -32°C	10°C~32°C	10°C~32°C
Shelves/Baskets	4/4	4/16	5/20	6/24	12/36	12/48	-
Drawers	-	-	-	--	-	-	12
Blood Storage Capacity	84	160	200	360	396	576	576
Door	Single-layer glass door				Double-layer glass door		
Casters	4	4	4	4	4	4	4
internal D (mm) (W*D*H)	455*370*1070	620*470*930	620*470*1200	680*620*1400	1030*590*1425	1320*700*1500	1320*700*1500
External Di(mm) (W*D*H)	560*630*1530	720*690*1530	720*690*1730	780*840*1960	1125*740*1980	1440*875*1980	1440*875*1980
Packing Di (mm)(W*D*H)	637*667*1647	797*787*1647	797*787*1857	857*937*2097	1254*840*2240	1600*1010*2100	1600*1010*2100
Net/gross Weight(kg)	107/120	146/160	158/174	204/227	227/279	310/360	310/360

Blood bags and Drawers can be added for all models



Model	Freeze 92	Freeze 188	Freeze 262	Freeze 428	Freeze 626	Freeze 100CH	Freeze 255CH	Freeze 380CH
Capacity(L)	92	188	262	428	626	100	255	380
Voltage/Frequency	220/50	220/50	220/50	220/50	220/50	220/50	220/50	220/50
Power(W)	220	290	340	1200	700	260	410	550
Temp. range(°C)	-10~-40	-10~-40	-10~-40	-10~-40	-20~-45	-20~-40	-20~-40	-20~-40
Drawer/Shelve	3/-	7/-	7/-	5/-	/3	/	/	/
Cabinet	Upright	Upright	Upright	Upright	Upright	Chest	Chest	Chest
Interior D (mm)(W*D*H)	465*440*645	380*420*1225	480*465*1430	760*630*1300	760*630*1310	500*385*625	1040*430*605	1376*457*625
Exterior D (mm)(W*D*H)	597*610*860	597*610*1475	657*685*1665	910*835*1990	1035*900*1980	680*600*805	1243*633*838	1550*633*838
Packing D (mm)(W*D*H)	659*699*934	659*699*1548	732*745*1756	1065*975*2125	1072*954*2130	735*645*855	1310*750*880	1630*755*880
Net/gross Weight(kg)	46/51	63/71	88/93	298/314	298/314	43/46	69/81	82/89

-40C° Deep Freezer



Model	UL288	UL388	UL486	UL490	UL628	UL728	UL826	ULT100	ULT420
Capacity(L)	288	388	486	490	628	728	826	100	420
Voltage/Frequency	220/50	220/50	220/50	220/50	220/50	220/50	220/50	220/50	220/50
Power(W)	930	1000	1000	1100	1200	1200	1200	650	1000
Temp. range(°C)	-40~-86	-40~-86	-40~-86	-40~-86	-40~-86	-40~-86	-40~-86	-40~-86	-40~-86
Ambient Temp.	10°C~32°C	10°C~32°C	10°C~32°C	10°C~32°C	10°C~32°C	10°C~32°C	10°C~32°C	10°C~32°C	10°C~32°C
Cabinet	Upright	Upright	Upright	Upright, double door	Upright	Upright	Upright	Chest	Chest
Shelves/Inner door	2/3	3/4	¾	3/4	3/4	3/4	3/4	-/1	-/3
Internal DI(mm)(W*D*H)	620*550*971	620*550*1310	590*630*1310	590*630*1310	760*630*1310	870*630*1310	870*710*1310	470*450*480	1367*462*652
Exterior DI(mm)(W*D*H)	915*870*1641	915*870*1980	945*900*1980	860*900*1980	1035*900*1980	1145*900*1980	1145*980*1980	769*825*1120	2130*870*1020
Packing DI (mm)(W*D*H)	952*924*1801	952*924*2130	982*954*2130	910*970*2130	1072*954*2130	1182*954*2130	1182*1034*2130	834*854*1235	2300*970*1164
Net/Gross Weight(kg)	217/236	255/276	290/310	295/335	301/323	350/380	380/410	138/160	310/357

-86°C Ultra Deep Freezer