



Model	Pharm260	Pharm360	Pharm610	Pharm890	Pharm1378
Capacity(L)	260	360	610	890	1378
Voltage/Frequency	220/50	220/50	220/50	220/50	220/50
Power(W)	340	430	490	900	800
Temp. range(°C)	2~8	2~8	2~8	2~8	2~8
Door	Single-layer glass	Single-layer glass	Single-layer glass	Double-layer glass	Double-layer glass
Shelves/Inner door	4/-	5/-	6/-	12/-	12/-
Internal DI(mm)(W*D*H)	550*460*1065	550*460*1340	680*620*1400	1030*590*1425	1320*700*1500
Exterior DI(mm)(W*D*H)	620*655*1720	620*655*1955	780*840*1960	1125*740*1980	1440*875*1980
Packing DI (mm)(W*D*H)	700*715*1810	700*715*2080	875*937*2097	1260*850*2150	1600*1010*2100
Net/Gross Weight(kg)	88/100	99/114	204/227	227/279	310/360

Model	Sun60
Storage Capacity(Litres)	21
Gross Capacity(Litres)	60
Ambient Temperature	10C~40C
Temp. Range (C)	2~8
Autonomy Time at 30C	130 h
Hold Over Time During Power Cut	110-120 h
Internal DI(mm)(W*D*H)	240*445*610
Exterior DI(mm)(W*D*H)	788*745*875
Net/Gross Weight(kg)	78/98