



Model	BB158	BB258	BB358	BB608	BB890	BB1308	BB1308K
Capacity(L)	158	258	358	608	890	1308	1308K
VAC /HZ	220/50	220/50	220/50	220/50	220/50	220/50	220/50
Power(W)	350	430	460	490	900	850	850
(°C)	4±1	4±1	4±1	4±1	4±1.5	4±1.5	4±1.5
Temp.control	Microprocessor control				Microprocessor control		
Ambient Temp.	10°C~38°C	10°C~38°C	10°C~38°C	10°C~38°C	10°C -32°C	10°C~32°C	10°C~32°C
Shelves/Baskets	4/4	4/16	5/20	6/24	12/36	12/48	-
Drawers	-	-	-	--	-	-	12
Blood Storage Capacity	84	160	200	360	396	576	576
Door	Single-layer glass door				Double-layer glass door		
Casters	4	4	4	4	4	4	4
internal D (mm) (W*D*H)	455*370*1070	620*470*930	620*470*1200	680*620*1400	1030*590*1425	1320*700*1500	1320*700*1500
External Di(mm) (W*D*H)	560*630*1530	720*690*1530	720*690*1730	780*840*1960	1125*740*1980	1440*875*1980	1440*875*1980
Packing Di (mm)(W*D*H)	637*667*1647	797*787*1647	797*787*1857	857*937*2097	1254*840*2240	1600*1010*2100	1600*1010*2100
Net/gross Weight(kg)	107/120	146/160	158/174	204/227	227/279	310/360	310/360

Blood bags and Drawers can be added for all models