

M930 Ion meter

pH Range: -2.000~20.000pH

pH Accuracy: ± 0.002 pH

Ion Range: 0.001~19999ppm

Ion Accuracy: $\pm 0.5\%$ F.S, $\pm 1\%$ F.S

Ion Concentration Range 0.001~19999ppm, mg/L, mol/L (Depending on measuring range of Ion Selective Electrode) Ion Concentration Accuracy $\pm 0.5\%$ Full Scale (Monovalent), $\pm 1\%$ Full Scale (Divalent)

Temperature Range: 0~105°C

Temperature Accuracy: $\pm 0.5^\circ\text{C}$

Temperature Compensation: 0~100°C

