



M900 Precision PH Meter

Model	M900-1	M900-2	M900-3
pH Range	0.00~14.00pH	0.00~14.00pH	0.00~14.00pH
pH Accuracy	±0.05pH	±0.01pH	±0.01pH
mV Range	-1999~1999mV	-1999~1999mV	-1999~1999mV
mV Accuracy	±10mV	±1mV	±1mV
Temperature Range	---	---	0~100°C
Temperature Accuracy	---	---	±1°C
Temperature Compensation	0~100°C, Manual	0~100°C, Manual	0~100°C, Automatic or Manual
Calibration Points	2 points		
Calibration Solutions	USA (pH4.01, 7.00, 10.01) or NIST (pH4.01, 6.86, 9.18)		
Power Requirements	DC9V, using AC adapters, 220VAC/50Hz		
Dimensions	210(L)×205(W)×65(H)mm		
Weight	1.7kg		